

VDS-150

Series

CroPharm


The Robust Entry into the world of Vertical Tablet Dedusting

For High Volume Pharmaceutical and
Nutraceutical Production



VDS-150

The Standard for
Vertical Tablet Dedusting



Combination Systems

Made in USA

Tablet Dedusting with In-line Metal Detection

Full mechanical separation of deduster and metal detector for highest detection accuracy

Ultra high-sensitivity metal detection with various reject options

Metal Detector controls
Integrated on one main frame

Customized solution to fit both deduster and metal detector on one rollaway stand



VDS-150 Technical Specifications

Specifications	VDS-150S	VDS-150M	VDS-150T
Inlet Height (in.)	30-36	30-36	30-36
Outlet Height (in.)	34-40	40-46	46-52
Avg. Output (tablets/hour)*	600,000*	600,000*	600,000*
Footprint (in.)	17 x 17	17 x 17	17 x 17
Product Contact Parts	SS316L	SS316L	SS316L
Outlet Positions	90°, 180°, 270°	90°, 180°, 270°	90°, 180°, 270°

VDS-150 Combination Technical Specifications (CEIA)

Specifications	VDS-150S-CEIA	VDS-150M-CEIA	VDS-150T-CEIA
Combination Unit			
Inlet Height (in.)	30-36	30-36	30-36
Outlet Height Combo (in.)	22-28	28-34	34-40
Avg. Output (tablets/hour)*	600,000*	600,000*	600,000*
Footprint (in.)	22 x 32	22 x 32	22 x 32
Detector Sensitivity	Ferrous	Non-Ferrous	Stainless Steel
CEIA THS/PH21E 100 x 40	0.3 mm	0.35 mm	0.5 mm
CEIA THS/PH21N 100 x 40	0.25 mm	0.3 mm	0.4 mm

Contact

USA (HQ)

Milford, CT USA

P: (203) 877-3859

E: sales@cropharm.com

W: cropharm.com

